
Revisions

SO 0191-5


Rev from	Rev to	ECO #
0608	2010-10	5253-10

Notes:


- BD Cat. Number 220245, 220246 & 220532
- Blank (Sheet) Size: Length: 8.5" Width: 11"
 Number of Pages: 56 Number of Sheets: 14
 Page Size: Length 8.5" Width 5.5" Final Folded Size: 8.5" x 5.5"
- Style (see illustrations below): # 5


#1


#2


#3


#4

#5

#6

#7

- See Specification Control Number OEM for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# 2755 Blue
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 7 Loveton Circle Sparks, MD 21152 USA	
Proofer	Date			
Checked By	Date			
Part Number: 54		Category and Description Package Insert, BD ESwab	Sheet: 1 of 57	A
			Scale: N/A	

BD Liquid Amies Elution Swab (ESwab) Collection and Transport System

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INTENDED USE

BD™ Liquid Amies Elution Swab (ESwab) Collection and Transport System is intended for the collection and transport of clinical specimens containing aerobes, anaerobes and fastidious bacteria from the collection site to the testing laboratory. In the laboratory, **BD** ESwab specimens are processed using standard clinical laboratory operating procedures for bacterial culture.

BD ESwab is supplied in a collection kit format. Each collection kit consists of a sterile package containing two components: a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 1 mL of Liquid Amies transport medium and a specimen collection swab which has a tip flocked with soft nylon fiber. Three collection kit configurations are available: one containing a regular size flocked nylon applicator; one containing a minitip size flocked nylon applicator; and one containing a flexible minitip size flocked nylon applicator.

SUMMARY AND EXPLANATION

One of the routine procedures in the diagnosis of bacteriological infections involves the collection and safe transportation of swab samples. This can be accomplished using the **BD** Liquid Amies Elution Swab (ESwab) Collection and Transport System. **BD** ESwab incorporates a modified Liquid Amies transporting medium, which can sustain the viability of a plurality of organisms that include clinically important aerobes, anaerobes and fastidious bacteria such as *Neisseria gonorrhoeae* during transit to the testing laboratory. The ESwab transport medium is a maintenance medium comprised of inorganic phosphate buffer, calcium and magnesium salts, and sodium chloride with a reduced environment due to the presence of sodium thioglycollate.¹

PRINCIPLES OF THE PROCEDURE

Once a swab sample is collected, it should be placed immediately into the **BD** ESwab transport tube where it comes into contact with the transport medium. Swab specimens for bacterial investigations collected using ESwab should be transported directly to the laboratory, preferably within 2 h of collection to maintain optimum organism viability.²⁻⁴ If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 h except for *Neisseria gonorrhoeae* cultures which should be processed within 24 h. Independent scientific studies on swab transport systems have shown that the viability of certain bacteria is superior at refrigerated temperatures compared with room temperature.⁵⁻⁹

REAGENTS

BD ESwab incorporates a modified Liquid Amies medium containing:

- Sodium chloride
- Potassium chloride
- Calcium chloride
- Magnesium chloride
- Monopotassium phosphate
- Disodium phosphate
- Sodium thioglycollate
- Distilled water

NOTE: The modified Liquid Amies Medium in **BD** ESwab transport tubes can have a cloudy appearance. This is normal and is due to the presence of salts in the medium formulation.

Warnings and Precautions

1. For *in vitro* Diagnostic Use.
2. Observe approved biohazard precautions and aseptic techniques. To be used only by adequately trained and qualified personnel.
3. Pathogenic microorganisms, including hepatitis viruses and Human Immunodeficiency Virus, may be present in clinical specimens. "Standard Precautions"¹⁰⁻¹³ and institutional guidelines should be followed in handling all items contaminated with blood and other body fluids.
4. Sterilize all biohazard waste including specimens, containers and media after their use.
5. Directions should be read and followed carefully.
6. Do not re-sterilize unused swabs.
7.  **BD** ESwab Collection and Transport System is for single use only; reuse may cause a risk of infection and/or inaccurate results.

8. Do not re-pack.
9. Not suitable to collect and transport microorganisms other than aerobes, anaerobes and fastidious bacteria.
10. Not suitable for any other application beyond the stated intended use.
11. The use of this product in association with a rapid diagnostic kit or with diagnostic instrumentation should be previously validated by the user.
12. Do not use if the swab is visibly damaged (i.e., if the swab tip or swab shaft is broken).
13. Do not use excessive force or pressure when collecting swab samples from patients as this may result in breakage of the swab shaft.
14. Do not ingest the medium.
15. Do not use the BD ESwab medium for pre-moistening or pre-wetting the applicator swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
16. BD ESwab should not be used if (1) there is evidence of damage or contamination to the product, (2) there is evidence of leakage, (3) the expiration date has passed, (4) the collection kit is open, or (5) there are other signs of deterioration.
17. Do not bend swab before use. If there is evidence to suggest the swab was bent do not use.

Storage and Handling

This product is ready for use and no further preparation is necessary. The product should be stored in its original packaging at 5 – 25°C until used. Do not incubate or freeze. Improper storage will result in a loss of efficacy. Do not use after expiration date, which is clearly printed on the outer box and on each individual sterile collection kit and the specimen transport tube label.

Materials Provided

BD ESwab is supplied in a collection kit format. Fifty (50) collection kits are contained in a shelf pack. Each collection kit consists of two components: a pre-labeled polypropylene screw-cap tube with conical shaped bottom filled with 1 mL of Liquid Amies transport medium, and a specimen collection swab which has a tip flocked with soft nylon fiber (see Figs 1 and 2). Three types of collection kit configurations are available: one containing a regular size flocked nylon swab applicator for sampling sites* (nose, throat, vagina and wounds); one containing a minitip size flocked nylon swab applicator for sampling sites* (eye, ear, nasal passages, nasopharynx, throat, urogenital tract and for pediatric sample collection); and one containing a flexible minitip size flocked nylon swab applicator for sampling sites* (nasopharynx and pediatric sample collection).

*Performance testing was not conducted using human specimens.

Fig 1. BD ESwab Collection Kit



Fig 2. BD ESwab Collection Kit Components

Materials Required But Not Provided

Appropriate materials for isolating and culturing aerobes, anaerobes and fastidious bacteria. These materials include culture media plates or tubes and incubation systems, gas jars or anaerobic workstations. Refer to laboratory reference manuals for recommended protocols for culture and identification techniques for aerobes, anaerobes and fastidious bacteria from clinical swab samples.¹⁴⁻²⁰

SPECIMEN COLLECTION, STORAGE AND TRANSPORT

Laboratory safety procedures should be followed. Caution should be used to minimize splashes and spills of medium during specimen collection and testing procedures described below.

Samples collected for bacteriological investigations which comprise the isolation of aerobes, anaerobes and fastidious bacteria such as *Nisseria gonorrhoeae* should be collected and handled following published manuals and guidelines.^{2,3,14-19}


To maintain optimum organism viability, transport specimens collected using BD ESwab directly to the laboratory, preferably within 2 h of collection.²⁻⁴ If immediate delivery or processing is delayed, then specimens should be refrigerated at 4 – 8°C or stored at room temperature (20 – 25°C) and processed within 48 h except for *Neisseria gonorrhoeae* cultures which should be processed within 24 h.

Specific requirements for the shipment and handling of specimens should be in full compliance with state and federal regulations.^{15,18,19} Shipping of specimens within medical institutions should comply with internal guidelines of the institution. All specimens should be processed as soon as they are received in the laboratory.

Specimen Collection

Proper specimen collection from the patient is extremely critical for successful isolation and identification of infectious organisms. For specific guidance regarding specimen collection procedures, consult published reference manuals.^{2,14,16-18,20}


NOTE: Swab should not be bent prior to sample collection.

1. Open the BD ESwab sample collection kit and remove the tube and the swab applicator.
2. Collect the specimen from the patient.
3. Aseptically unscrew and remove the cap from the tube.
4. Insert the swab into the tube and bend the swab shaft at the breakpoint indicated by the colored line marked on the swab shaft against the tube to break the shaft. Discard the broken handle part of the swab shaft into an approved medical waste disposal container.
5. Replace cap on the tube and secure tightly.
6. Write patient information on the tube label or apply patient identification label. Send the specimen to the test laboratory.


During specimen collection when handling the swab applicator, the operator must not touch the area below the marked breakpoint indication line; that is the area from the line to the tip of the nylon flocked swab (see Fig 3), as this will lead to contamination of the applicator shaft and the culture thus invalidating the test results.

Fig 3. Collection swab showing breakpoint indication line and area for holding the applicator

After the swab sample is taken from the patient, the swab applicator shaft is broken off at the colored breakpoint indication line by bending the shaft on the side of the tube allowing the swab to fall into the BD ESwab tube of transport medium. The operator then discards the handle part of the swab into an approved medical waste disposal container. The tube's screw cap is then replaced and secured tightly. The action of screwing the cap onto the tube moves the end of the broken swab shaft into a funnel shaped molded docking receptacle in the cap (see Fig 4). This molded funnel shape captures the end of the broken applicator shaft and secures it firmly in the dock by friction grip.

Fig 4. Capture of broken swab applicator stick by BD ESwab tube cap

In the testing laboratory when the BD ESwab cap is unscrewed and removed, the swab applicator stick is securely attached to the cap. This feature allows the operator to conveniently remove the swab and perform various microbiology analyses using the tube cap as a handle to hold and manipulate the swab. Due to the flexibility of the shaft of the minitip swabs (220246 & 220532), the capture cap feature is not applicable, as the broken applicator may not firmly fit in the cap. Use sterile forceps to extract the applicator from the tube or from the cap if the swab has been partially captured.

TEST PROCEDURE

Plating BD ESwab Specimen Cultures in the Laboratory

BD ESwab specimens should be processed for bacteriological culture using recommended culture media and laboratory techniques which will depend on the specimen type and the organism under investigation. For recommended culture media and techniques for the isolation and identification of bacteria from clinical swab specimens refer to published microbiology manuals and guidelines.^{14,17,20-22}

Culture investigations of swab specimens for the presence of aerobic bacteria, anaerobic bacteria and fastidious bacteria such as *Neisseria gonorrhoeae* routinely involve the use of solid agar culture medium in Petri dish plates. The procedure for inoculation of BD ESwab samples onto solid agar in Petri dishes is as follows.

1. Vigorously shake the BD ESwab tube containing the swab sample between the thumb and forefinger for 5 s or mix the tube using a vortex mixer for 5 s to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the liquid transport medium.
2. Unscrew the BD ESwab cap to remove the swab applicator.

3. Roll the tip of the **BD ESwab** applicator onto the surface of one quadrant of the culture media plate to provide the primary inoculum.
4. If it is necessary to culture the swab specimen onto additional culture media plates, return the **BD ESwab** applicator to the transport medium tube for two seconds to absorb and recharge the applicator tip with transport medium/patient sample suspension then repeat Step No. 3.


The procedure described above utilizes the **BD ESwab** applicator like an inoculation wand to transfer the suspension of patient sample in transport medium onto the surface of a culture plate creating the primary inoculum (see Fig 5). Alternatively, the operator can vortex or mix the **BD ESwab** tube with the swab inside for 5 s and then transfer 100 μ L volumes of the suspension onto each culture plate using a volumetric pipetor and sterile pipet tips. Standard laboratory techniques should then be used to streak the primary inoculum of patient sample across the surface of the culture plate (see Fig 6).

Fig 5. Procedures for inoculation of BD ESwab specimens onto solid agar in Petri dishes

1. Using swab to inoculate specimen
2. Using pipetor and sterile pipet tips to inoculate 100 μ L of specimen

Fig 6. Procedure for streaking BD ESwab specimens on agar Petri dishes for primary isolation²³

Seed a primary inoculum of **BD ESwab** specimen onto the surface of an appropriate culture media following recommended microbiology procedures.

Use a sterile bacteriology loop to streak the primary inoculum across the surface to isolate organisms on the culture media.

Preparation of Gram Stain Smears of BD ESwab Specimens

Laboratory analysis of clinical swab samples collected from certain sites on the patient can routinely include microscopic examination of stained preparations ("direct smears") using the Gram stain procedure. This can provide valuable information to physicians who are managing patients with infectious diseases.²⁴ There are many instances in which a Gram stain can assist in making a diagnosis; for example, with swabs taken from the endocervix or male urethra to investigate suspected *Neisseria gonorrhoeae* infections or vaginal swabs to diagnose bacterial vaginosis.²⁵⁻³⁰ The Gram stain can also help to judge specimen quality and contribute to the selection of culture media especially with mixed flora.³¹

Microscope slides of patient specimens transported in **BD ESwab** transport system can be prepared for Gram stain analysis, as describe below, by sampling an aliquot of vortexed suspension of the swab.^{17,31}

1. Take a clean glass microscope slide, place it on a flat surface and inscribe an area using a diamond-tipped or similar glass marker to identify the location of the specimen inoculum.
2. Vortex or mix the **BD ESwab** tube containing the swab sample for 5 s to release the sample from the swab tip and evenly disperse and suspend the patient specimen in the Liquid Amies transport medium.
3. Unscrew the **BD ESwab** cap and using a sterile pipet, transfer 1 – 2 drops of Liquid Amies medium to the inscribed area on the glass slide.
4. Allow the specimen on the slide to air dry or place the slide in a 60°C electric slide warmer.
5. Smears may be fixed with heat or methanol. Methanol fixation may be preferable as it prevents lysis of red blood cells, avoids damage to all host cells and results in a cleaner background.^{17,24,31}
6. Follow published laboratory reference manuals and guidelines for performing the Gram stain. If commercial Gram stain reagents are used, it is important to comply with instructions in the manufacturer's product insert for performance test procedure.

For further information or guidance on the preparation of specimen slides for microscopic analysis, for information on Gram staining procedures and the interpretation and reporting of microscopic analysis, consult published laboratory reference manuals.^{16,21,22,24,31}

QUALITY CONTROL

All lot numbers of the **BD ESwab** are tested for sterility and swab applicators are tested to ensure they are non-toxic to bacteria. **BD ESwab** Liquid Amies transport medium is tested for pH stability and bio-burden using Gram stain microscopic examination to ensure acceptable levels as defined in the Clinical and Laboratory Standards Institute (CLSI) Approved Standard M40-A.⁴ Each production lot of **BD ESwab** is quality control tested before release for ability to maintain viable bacteria at both refrigerated temperatures (4 – 8°C) and room temperature (20 – 25°C) for specified time points with a panel of aerobes, anaerobes and fastidious bacteria using both Roll-Plate and Swab Elution Methods.⁴ Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C) which should correspond to ≤ 1 log increase in growth at a specified time point. Procedures for quality control of bacteriology transport devices using a quantitative Swab Elution Method and qualitative Roll-Plate

Method are described in CLSI M40-A and other publications.^{4,5,7,8,32-34} If aberrant quality control results are noted, patient results should not be reported.

LIMITATIONS OF THE PROCEDURE

1. Condition, timing, and volume of specimen collected for culture are significant variables in obtaining reliable culture results. Follow recommended guidelines for specimen collection.^{2,3,14,16,17,20,21}
2. **BD** ESwab is intended for use as a collection and transport medium for aerobes, anaerobes and fastidious bacteria such as *Neisseria gonorrhoeae*.
3. **BD** ESwab Collection and Transport System is intended to be used with the medium tubes and swabs provided in the kit. The use of tubes of medium or swabs from any other source are not qualified for use with **BD** ESwab and could affect the performance of the product and laboratory test results.

EXPECTED VALUES

Results obtained will largely depend on proper and adequate specimen collection, as well as timely transport and processing in the laboratory.

PERFORMANCE CHARACTERISTICS

In the routine clinical laboratory, the Roll-Plate Method is the primary means of inoculating swab transport devices onto plated media. A limitation of the Roll-Plate Method⁴ for bacterial viability performance testing is that it is not a quantitative method; it is, at best, a semiquantitative approximation. On the other hand, quantitative viability performance methods such as the Swab Elution Method⁴ do not reflect the standard protocol used in most clinical laboratories. Whereas the Swab Elution Method allows a quantitative measurement of the ability of a transport system to maintain viable organisms, the Roll-Plate technique takes into consideration some mechanical variables of the direct swabbing action that exist in the clinical laboratory, and which can influence the release of the sample onto culture plates. Because of this, both methods of performing viability studies were used to determine the performance characteristics of the **BD** ESwab Collection and Transport System.

The test procedures employed for determining bacterial viability performance were based upon the quality control methods described in CLSI M40-A.^{4,5,7,8,32-34} The test organisms utilized in this study were those specifically prescribed in CLSI M40-A for establishing performance claims and quality control of swab transport systems and include a representative panel of aerobes, anaerobes and fastidious bacteria. An additional group of organisms not required or specified by CLSI M40-A were tested in order to provide further information on the survival of specific bacteria. Bacterial viability studies were performed on the **BD** ESwab at two different ranges of temperature, 4 – 8 °C and 20 – 25 °C, corresponding to refrigerator and room temperature, respectively. Swabs accompanying each transport system were inoculated in triplicate with 100 µL of specific concentrations of organism suspension. Swabs were then placed in their respective transport medium tubes and were held for 0 h, 24 h and 48 h. At the appropriate time intervals, each swab was processed according to the Roll-Plate or Swab Elution Method.

Organisms evaluated were divided into three main groups (see note below):

1. Aerobes and Facultative Anaerobes:

Pseudomonas aeruginosa ATCC™ BAA-427, *Streptococcus pyogenes* ATCC 19615, *Streptococcus pneumoniae* ATCC 6305, *Haemophilus influenzae* ATCC 10211.

2. Anaerobes:

Bacteroides fragilis ATCC 25285, *Peptostreptococcus anaerobius* ATCC 27337, *Fusobacterium nucleatum* ATCC 25586, *Propionibacterium acnes* ATCC 6919, *Prevotella melaninogenica* ATCC 25845.

3. Fastidious Bacteria:

Neisseria gonorrhoeae ATCC 43069.

Additional organisms evaluated:

Enterococcus faecalis (vancomycin-resistant *Enterococcus*, VRE) ATCC 51299, *Staphylococcus aureus* (methicillin-resistant *Staphylococcus aureus*, MRSA) ATCC 43300, *Streptococcus agalactiae* (Group B *Streptococcus*) ATCC 13813, *Clostridium perfringens* ATCC 13124, *Clostridium sporogenes* ATCC 3584, *Fusobacterium necrophorum* ATCC 25286, *Peptococcus magnus* ATCC 29328.

The results for the bacterial strains tested using the BD ESwab System are shown in the following tables.

**SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES
ROLL-PLATE METHOD, 4-8°C**

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	BD ESwab Lot	Average of CFUs recovered at 0 h	Average of CFUs recovered at 24 h	Average of CFUs recovered at 48 h	Interpretation
<i>Pseudomonas aeruginosa</i> ATCC BAA-427	diluted 10 ^{-3.5}	1	261.7	210.7	59.3	Acceptable Recovery
		2	258.3	206.3	54.7	Acceptable Recovery
		3	268.0	203.3	56.7	Acceptable Recovery
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 10 ⁻³	1	292.7	142.0	49.0	Acceptable Recovery
		2	283.6	138.3	49.3	Acceptable Recovery
		3	285.6	145.3	48.0	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 10 ^{-1.5}	1	193.3	60.7	29.7	Acceptable Recovery
		2	194.7	61.7	32.3	Acceptable Recovery
		3	196.7	64.0	35.0	Acceptable Recovery
<i>Haemophilus influenzae</i> ATCC 10211	diluted 10 ^{-3.5}	1	277.7	121.0	27.3	Acceptable Recovery
		2	267.7	111.3	19.7	Acceptable Recovery
		3	260.7	101.3	17.3	Acceptable Recovery
<i>Bacteroides fragilis</i> ATCC 25285	diluted 10 ⁻³	1	288.3	93.7	54.0	Acceptable Recovery
		2	278.3	83.7	44.0	Acceptable Recovery
		3	272.7	74.3	29.7	Acceptable Recovery
<i>Peptostreptococcus anaerobius</i> ATCC 27337	diluted 10 ^{-2.5}	1	286.7	180.3	22.7	Acceptable Recovery
		2	290.0	182.7	21.3	Acceptable Recovery
		3	284.3	187.3	23.3	Acceptable Recovery
<i>Fusobacterium nucleatum</i> ATCC 25586	diluted 10 ^{-1.5}	1	272.0	110.0	19.0	Acceptable Recovery
		2	275.0	102.0	16.7	Acceptable Recovery
		3	272.0	111.0	22.0	Acceptable Recovery
<i>Propionibacterium acnes</i> ATCC 6919	diluted 10 ⁻³	1	290.7	156.7	48.7	Acceptable Recovery
		2	288.3	151.3	40.7	Acceptable Recovery
		3	290.7	154.7	47.0	Acceptable Recovery
<i>Prevotella melaninogenica</i> ATCC 25845	diluted 10 ^{-2.5}	1	292.3	169.3	29.3	Acceptable Recovery
		2	288.0	168.3	31.0	Acceptable Recovery
		3	292.7	169.7	29.7	Acceptable Recovery
<i>Neisseria gonorrhoeae</i> ATCC 43069	diluted 10 ⁻³	1	234.7	19.7		Acceptable Recovery
		2	244.7	24.3		Acceptable Recovery
		3	246.3	23.7		Acceptable Recovery
<i>Enterococcus faecalis</i> (VRE) ATCC 51299	diluted 10 ^{-3.5}	1	240.0	109.3	41.3	Acceptable Recovery
		2	230.0	101.7	37.3	Acceptable Recovery
		3	247.7	102.3	41.0	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted 10 ^{-3.5}	1	238.0	98.0	50.3	Acceptable Recovery
		2	238.7	98.7	49.0	Acceptable Recovery
		3	236.3	96.3	48.0	Acceptable Recovery
<i>Streptococcus agalactiae</i> (Group B Strep) ATCC 13813	diluted 10 ^{-3.5}	1	290.0	116.7	56.3	Acceptable Recovery
		2	292.3	116.7	58.3	Acceptable Recovery
		3	291.0	116.3	56.7	Acceptable Recovery
<i>Clostridium perfringens</i> ATCC 13124	diluted 10 ^{-3.5}	1	283.3	162.0	48.7	Acceptable Recovery
		2	279.3	152.0	41.7	Acceptable Recovery
		3	273.3	145.3	44.0	Acceptable Recovery
<i>Clostridium sporogenes</i> ATCC 3584	diluted 10 ^{-3.5}	1	248.3	100.3	43.7	Acceptable Recovery
		2	247.0	94.7	38.3	Acceptable Recovery
		3	238.3	91.3	33.7	Acceptable Recovery
<i>Fusobacterium necrophorum</i> ATCC 25286	diluted 10 ^{-2.5}	1	288.0	146.7	51.3	Acceptable Recovery
		2	278.0	136.7	41.3	Acceptable Recovery
		3	274.7	132.7	47.7	Acceptable Recovery
<i>Peptococcus magnus</i> ATCC 29328	diluted 10 ^{-2.5}	1	284.3	153.7	42.3	Acceptable Recovery
		2	288.0	152.3	43.3	Acceptable Recovery
		3	274.3	144.3	34.0	Acceptable Recovery

SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES
ROLL-PLATE METHOD, 20-25°C

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	BD ESwab Lot	Average of CFUs recovered at 0 h	Average of CFUs recovered at 24 h	Average of CFUs recovered at 48 h	Interpretation
<i>Pseudomonas aeruginosa</i> ATCC BAA-427	diluted 10 ^{-3.5}	1	261.7	190.0	51.7	Acceptable Recovery
		2	258.3	178.0	44.7	Acceptable Recovery
		3	268.0	192.3	49.0	Acceptable Recovery
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 10 ⁻³	1	292.7	108.0	33.0	Acceptable Recovery
		2	283.6	115.7	33.0	Acceptable Recovery
		3	285.6	109.7	31.0	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 10 ^{-1.5}	1	193.3	56.0	23.0	Acceptable Recovery
		2	194.7	54.7	21.7	Acceptable Recovery
		3	196.7	58.7	22.0	Acceptable Recovery
<i>Haemophilus influenzae</i> ATCC 10211	diluted 10 ^{-3.5}	1	277.7	113.3	19.3	Acceptable Recovery
		2	267.7	98.3	17.0	Acceptable Recovery
		3	260.7	88.3	11.0	Acceptable Recovery
<i>Bacteroides fragilis</i> ATCC 25285	diluted 10 ⁻³	1	288.3	76.3	40.7	Acceptable Recovery
		2	278.3	67.7	32.7	Acceptable Recovery
		3	272.7	60.7	26.7	Acceptable Recovery
<i>Peptostreptococcus anaerobius</i> ATCC 27337	diluted 10 ^{-2.5}	1	286.7	164.0	14.3	Acceptable Recovery
		2	290.0	154.0	14.0	Acceptable Recovery
		3	284.3	164.0	15.7	Acceptable Recovery
<i>Fusobacterium nucleatum</i> ATCC 25586	diluted 10 ^{-1.5}	1	272.0	86.3	17.3	Acceptable Recovery
		2	275.0	78.0	12.7	Acceptable Recovery
		3	272.0	76.3	17.3	Acceptable Recovery
<i>Propionibacterium acnes</i> ATCC 6919	diluted 10 ⁻³	1	290.7	107.3	36.0	Acceptable Recovery
		2	288.3	97.3	28.3	Acceptable Recovery
		3	290.7	105.3	34.7	Acceptable Recovery
<i>Prevotella melaninogenica</i> ATCC 25845	diluted 10 ^{-2.5}	1	292.3	92.3	16.7	Acceptable Recovery
		2	288.0	93.3	15.0	Acceptable Recovery
		3	292.7	92.3	17.3	Acceptable Recovery
<i>Neisseria gonorrhoeae</i> ATCC 43069	diluted 10 ⁻³	1	234.7	13.7		Acceptable Recovery
		2	244.7	15.7		Acceptable Recovery
		3	246.3	18.0		Acceptable Recovery
<i>Enterococcus faecalis</i> (VRE) ATCC 51299	diluted 10 ^{-3.5}	1	240.0	93.7	32.7	Acceptable Recovery
		2	230.0	89.0	27.7	Acceptable Recovery
		3	247.7	86.0	29.3	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted 10 ^{-3.5}	1	238.0	74.3	44.0	Acceptable Recovery
		2	238.7	73.3	42.7	Acceptable Recovery
		3	236.3	76.3	42.3	Acceptable Recovery
<i>Streptococcus agalactiae</i> (Group B Strep) ATCC 13813	diluted 10 ^{-3.5}	1	290.0	88.0	47.7	Acceptable Recovery
		2	292.3	87.0	46.0	Acceptable Recovery
		3	291.0	86.3	46.3	Acceptable Recovery
<i>Clostridium perfringens</i> ATCC 13124	diluted 10 ^{-3.5}	1	283.3	110.7	37.0	Acceptable Recovery
		2	279.3	99.7	32.0	Acceptable Recovery
		3	273.3	92.0	32.0	Acceptable Recovery
<i>Clostridium sporogenes</i> ATCC 3584	diluted 10 ^{-3.5}	1	248.3	91.3	36.0	Acceptable Recovery
		2	247.0	86.3	31.7	Acceptable Recovery
		3	238.3	73.3	29.0	Acceptable Recovery
<i>Fusobacterium necrophorum</i> ATCC 25286	diluted 10 ^{-2.5}	1	288.0	107.3	40.3	Acceptable Recovery
		2	278.0	97.3	30.3	Acceptable Recovery
		3	274.7	97.0	33.7	Acceptable Recovery
<i>Peptococcus magnus</i> ATCC 29328	diluted 10 ^{-2.5}	1	284.3	107.3	31.3	Acceptable Recovery
		2	288.0	106.7	31.0	Acceptable Recovery
		3	274.3	97.3	24.3	Acceptable Recovery

SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES
SWAB ELUTION METHOD, 4-8°C

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	BD ESwab Lot	Average of CFUs recovered at 0 h	Average of CFUs recovered at 24 h	Average of CFUs recovered at 48 h	Log ₁₀ decline	Interpretation
<i>Pseudomonas aeruginosa</i> ATCC BAA-427	diluted 1:10	1	1.4 x 10 ⁶	1.1 x 10 ⁶	2.7 x 10 ⁵	-0.71	Acceptable Recovery
		2	1.4 x 10 ⁶	1.0 x 10 ⁶	2.6 x 10 ⁵	-0.73	Acceptable Recovery
		3	1.5 x 10 ⁶	9.7 x 10 ⁵	2.6 x 10 ⁵	-0.76	Acceptable Recovery
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 1:10	1	6.0 x 10 ⁵	2.9 x 10 ⁵	6.0 x 10 ⁴	-1.00	Acceptable Recovery
		2	6.0 x 10 ⁵	2.9 x 10 ⁵	6.5 x 10 ⁴	-0.97	Acceptable Recovery
		3	6.1 x 10 ⁵	3.0 x 10 ⁵	6.8 x 10 ⁴	-0.95	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 1:10	1	1.8 x 10 ⁶	6.0 x 10 ⁵	2.0 x 10 ⁵	-0.95	Acceptable Recovery
		2	1.8 x 10 ⁶	6.9 x 10 ⁵	2.0 x 10 ⁵	-0.95	Acceptable Recovery
		3	1.8 x 10 ⁶	6.4 x 10 ⁵	1.9 x 10 ⁵	-0.98	Acceptable Recovery
<i>Haemophilus influenzae</i> ATCC 10211	diluted 1:10	1	3.9 x 10 ⁶	9.6 x 10 ⁵	3.9 x 10 ⁵	-1.00	Acceptable Recovery
		2	3.8 x 10 ⁶	9.9 x 10 ⁵	3.6 x 10 ⁵	-1.02	Acceptable Recovery
		3	3.7 x 10 ⁶	8.9 x 10 ⁵	2.8 x 10 ⁵	-1.12	Acceptable Recovery
<i>Bacteroides fragilis</i> ATCC 25285	diluted 1:10	1	8.6 x 10 ⁵	3.7 x 10 ⁵	1.5 x 10 ⁵	-0.76	Acceptable Recovery
		2	8.4 x 10 ⁵	3.5 x 10 ⁵	1.4 x 10 ⁵	-0.78	Acceptable Recovery
		3	8.2 x 10 ⁵	3.3 x 10 ⁵	1.3 x 10 ⁵	-0.80	Acceptable Recovery
<i>Peptostreptococcus anaerobius</i> ATCC 27337	diluted 1:10	1	1.6 x 10 ⁶	9.7 x 10 ⁵	1.2 x 10 ⁵	-1.12	Acceptable Recovery
		2	1.7 x 10 ⁶	9.6 x 10 ⁵	1.1 x 10 ⁵	-1.16	Acceptable Recovery
		3	1.7 x 10 ⁶	9.5 x 10 ⁵	1.1 x 10 ⁵	-1.19	Acceptable Recovery
<i>Fusobacterium nucleatum</i> ATCC 25586	diluted 1:10	1	2.4 x 10 ⁶	7.0 x 10 ⁵	1.8 x 10 ⁵	-1.12	Acceptable Recovery
		2	2.4 x 10 ⁶	6.9 x 10 ⁵	1.8 x 10 ⁵	-1.12	Acceptable Recovery
		3	2.4 x 10 ⁶	6.8 x 10 ⁵	1.9 x 10 ⁵	-1.10	Acceptable Recovery
<i>Propionibacterium acnes</i> ATCC 6919	diluted 1:10	1	3.8 x 10 ⁶	1.9 x 10 ⁶	6.9 x 10 ⁵	-0.74	Acceptable Recovery
		2	3.7 x 10 ⁶	1.8 x 10 ⁶	6.0 x 10 ⁵	-0.79	Acceptable Recovery
		3	3.7 x 10 ⁶	1.8 x 10 ⁶	5.9 x 10 ⁵	-0.80	Acceptable Recovery
<i>Prevotella melaninogenica</i> ATCC 25845	diluted 1:10	1	3.1 x 10 ⁶	9.3 x 10 ⁵	2.7 x 10 ⁵	-1.06	Acceptable Recovery
		2	3.0 x 10 ⁶	9.3 x 10 ⁵	2.7 x 10 ⁵	-1.05	Acceptable Recovery
		3	3.2 x 10 ⁶	9.3 x 10 ⁵	2.6 x 10 ⁵	-1.09	Acceptable Recovery
<i>Neisseria gonorrhoeae</i> ATCC 43069	diluted 1:10	1	3.6 x 10 ⁶	2.8 x 10 ⁵		-1.11	Acceptable Recovery
		2	3.5 x 10 ⁶	2.7 x 10 ⁵		-1.11	Acceptable Recovery
		3	3.4 x 10 ⁶	2.5 x 10 ⁵		-1.13	Acceptable Recovery
<i>Enterococcus faecalis</i> (VRE) ATCC 51299	diluted 1:10	1	1.4 x 10 ⁶	8.4 x 10 ⁵	2.5 x 10 ⁵	-0.75	Acceptable Recovery
		2	1.4 x 10 ⁶	8.2 x 10 ⁵	2.5 x 10 ⁵	-0.75	Acceptable Recovery
		3	1.4 x 10 ⁶	8.5 x 10 ⁵	2.6 x 10 ⁵	-0.73	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted 1:10	1	9.9 x 10 ⁵	7.7 x 10 ⁵	1.9 x 10 ⁵	-0.72	Acceptable Recovery
		2	9.8 x 10 ⁵	7.6 x 10 ⁵	1.8 x 10 ⁵	-0.73	Acceptable Recovery
		3	1.0 x 10 ⁶	7.6 x 10 ⁵	2.0 x 10 ⁵	-0.70	Acceptable Recovery
<i>Streptococcus agalactiae</i> (Group B Strep) ATCC 13813	diluted 1:10	1	5.5 x 10 ⁶	3.4 x 10 ⁶	8.1 x 10 ⁵	-0.83	Acceptable Recovery
		2	5.6 x 10 ⁶	3.6 x 10 ⁶	8.0 x 10 ⁵	-0.85	Acceptable Recovery
		3	5.4 x 10 ⁶	3.4 x 10 ⁶	7.8 x 10 ⁵	-0.84	Acceptable Recovery
<i>Clostridium perfringens</i> ATCC 13124	diluted 1:10	1	2.3 x 10 ⁶	1.3 x 10 ⁶	3.9 x 10 ⁵	-0.77	Acceptable Recovery
		2	2.3 x 10 ⁶	1.2 x 10 ⁶	3.6 x 10 ⁵	-0.81	Acceptable Recovery
		3	2.2 x 10 ⁶	1.2 x 10 ⁶	3.2 x 10 ⁵	-0.84	Acceptable Recovery
<i>Clostridium sporogenes</i> ATCC 3584	diluted 1:10	1	6.5 x 10 ⁵	3.0 x 10 ⁵	1.2 x 10 ⁵	-0.73	Acceptable Recovery
		2	6.4 x 10 ⁵	4.0 x 10 ⁵	1.2 x 10 ⁵	-0.73	Acceptable Recovery
		3	6.4 x 10 ⁵	2.9 x 10 ⁵	1.1 x 10 ⁵	-0.76	Acceptable Recovery
<i>Fusobacterium necrophorum</i> ATCC 25286	diluted 1:10	1	9.6 x 10 ⁵	4.2 x 10 ⁵	1.7 x 10 ⁵	-0.75	Acceptable Recovery
		2	9.7 x 10 ⁵	4.3 x 10 ⁵	1.8 x 10 ⁵	-0.73	Acceptable Recovery
		3	9.4 x 10 ⁵	4.1 x 10 ⁵	1.6 x 10 ⁵	-0.77	Acceptable Recovery
<i>Peptococcus magnus</i> ATCC 29328	diluted 1:10	1	4.9 x 10 ⁶	2.9 x 10 ⁶	8.6 x 10 ⁵	-0.76	Acceptable Recovery
		2	4.9 x 10 ⁶	2.8 x 10 ⁶	8.7 x 10 ⁵	-0.75	Acceptable Recovery
		3	4.8 x 10 ⁶	2.8 x 10 ⁶	7.9 x 10 ⁵	-0.78	Acceptable Recovery

SUMMARY OF RESULTS FOR BACTERIAL RECOVERY STUDIES
SWAB ELUTION METHOD, 20-25°C

Organism	Dilution: 0.5 McFarland bacterial suspension with saline	BD ESwab Lot	Average of CFUs recovered at 0 h	Average of CFUs recovered at 24 h	Average of CFUs recovered at 48 h	Log ₁₀ decline	Interpretation
<i>Pseudomonas aeruginosa</i> ATCC BAA-427	diluted 1:10	1	1.4 x 10 ⁶	9.8 x 10 ⁵	2.7 x 10 ⁵	-0.71	Acceptable Recovery
		2	1.4 x 10 ⁶	9.6 x 10 ⁵	2.5 x 10 ⁵	-0.75	Acceptable Recovery
		3	1.5 x 10 ⁶	9.8 x 10 ⁵	2.3 x 10 ⁵	-0.81	Acceptable Recovery
<i>Streptococcus pyogenes</i> ATCC 19615	diluted 1:10	1	6.0 x 10 ⁵	2.6 x 10 ⁵	4.5 x 10 ⁴	-1.12	Acceptable Recovery
		2	6.0 x 10 ⁵	2.5 x 10 ⁵	4.1 x 10 ⁴	-1.17	Acceptable Recovery
		3	6.1 x 10 ⁵	2.5 x 10 ⁵	4.2 x 10 ⁴	-1.16	Acceptable Recovery
<i>Streptococcus pneumoniae</i> ATCC 6305	diluted 1:10	1	1.8 x 10 ⁶	4.4 x 10 ⁵	1.6 x 10 ⁵	-1.05	Acceptable Recovery
		2	1.8 x 10 ⁶	4.7 x 10 ⁵	1.5 x 10 ⁵	-1.08	Acceptable Recovery
		3	1.8 x 10 ⁶	4.7 x 10 ⁵	1.5 x 10 ⁵	-1.08	Acceptable Recovery
<i>Haemophilus influenzae</i> ATCC 10211	diluted 1:10	1	3.9 x 10 ⁶	8.2 x 10 ⁵	3.2 x 10 ⁵	-1.09	Acceptable Recovery
		2	3.8 x 10 ⁶	8.2 x 10 ⁵	2.9 x 10 ⁵	-1.12	Acceptable Recovery
		3	3.7 x 10 ⁶	7.2 x 10 ⁵	2.2 x 10 ⁵	-1.23	Acceptable Recovery
<i>Bacteroides fragilis</i> ATCC 25285	diluted 1:10	1	8.6 x 10 ⁵	3.8 x 10 ⁵	1.2 x 10 ⁵	-0.86	Acceptable Recovery
		2	8.4 x 10 ⁵	3.7 x 10 ⁵	1.2 x 10 ⁵	-0.85	Acceptable Recovery
		3	8.2 x 10 ⁵	3.5 x 10 ⁵	1.0 x 10 ⁵	-0.91	Acceptable Recovery
<i>Peptostreptococcus anaerobius</i> ATCC 27337	diluted 1:10	1	1.6 x 10 ⁶	8.5 x 10 ⁵	1.1 x 10 ⁵	-1.16	Acceptable Recovery
		2	1.7 x 10 ⁶	8.5 x 10 ⁵	9.9 x 10 ⁴	-1.23	Acceptable Recovery
		3	1.7 x 10 ⁶	8.3 x 10 ⁵	9.8 x 10 ⁴	-1.24	Acceptable Recovery
<i>Fusobacterium nucleatum</i> ATCC 25586	diluted 1:10	1	2.4 x 10 ⁶	6.6 x 10 ⁵	1.6 x 10 ⁵	-1.18	Acceptable Recovery
		2	2.4 x 10 ⁶	6.4 x 10 ⁵	1.6 x 10 ⁵	-1.18	Acceptable Recovery
		3	2.4 x 10 ⁶	6.5 x 10 ⁵	1.7 x 10 ⁵	-1.15	Acceptable Recovery
<i>Propionibacterium acnes</i> ATCC 6919	diluted 1:10	1	3.8 x 10 ⁶	1.3 x 10 ⁶	4.3 x 10 ⁵	-0.95	Acceptable Recovery
		2	3.7 x 10 ⁶	1.2 x 10 ⁶	3.3 x 10 ⁵	-1.05	Acceptable Recovery
		3	3.7 x 10 ⁶	1.2 x 10 ⁶	3.4 x 10 ⁵	-1.04	Acceptable Recovery
<i>Prevotella melaninogenica</i> ATCC 25845	diluted 1:10	1	3.1 x 10 ⁶	5.9 x 10 ⁵	2.1 x 10 ⁵	-1.17	Acceptable Recovery
		2	3.0 x 10 ⁶	5.9 x 10 ⁵	2.1 x 10 ⁵	-1.15	Acceptable Recovery
		3	3.2 x 10 ⁶	6.0 x 10 ⁵	2.1 x 10 ⁵	-1.18	Acceptable Recovery
<i>Neisseria gonorrhoeae</i> ATCC 43069	diluted 1:10	1	3.6 x 10 ⁶	2.2 x 10 ⁵		-1.21	Acceptable Recovery
		2	3.5 x 10 ⁶	2.1 x 10 ⁵		-1.22	Acceptable Recovery
		3	3.4 x 10 ⁶	1.9 x 10 ⁵		-1.25	Acceptable Recovery
<i>Enterococcus faecalis</i> (VRE) ATCC 51299	diluted 1:10	1	1.4 x 10 ⁶	7.6 x 10 ⁵	2.1 x 10 ⁵	-0.82	Acceptable Recovery
		2	1.4 x 10 ⁶	7.5 x 10 ⁵	2.0 x 10 ⁵	-0.85	Acceptable Recovery
		3	1.4 x 10 ⁶	7.5 x 10 ⁵	1.9 x 10 ⁵	-0.87	Acceptable Recovery
<i>Staphylococcus aureus</i> (MRSA) ATCC 43300	diluted 1:10	1	9.9 x 10 ⁵	6.9 x 10 ⁵	1.1 x 10 ⁵	-0.95	Acceptable Recovery
		2	9.8 x 10 ⁵	6.5 x 10 ⁵	1.2 x 10 ⁵	-0.91	Acceptable Recovery
		3	1.0 x 10 ⁶	6.6 x 10 ⁵	1.2 x 10 ⁵	-0.92	Acceptable Recovery
<i>Streptococcus agalactiae</i> (Group B Strep) ATCC 13813	diluted 1:10	1	5.5 x 10 ⁶	3.4 x 10 ⁶	5.4 x 10 ⁵	-1.01	Acceptable Recovery
		2	5.6 x 10 ⁶	3.3 x 10 ⁶	5.4 x 10 ⁵	-1.02	Acceptable Recovery
		3	5.4 x 10 ⁶	3.6 x 10 ⁶	5.5 x 10 ⁵	-0.99	Acceptable Recovery
<i>Clostridium perfringens</i> ATCC 13124	diluted 1:10	1	2.3 x 10 ⁶	1.0 x 10 ⁶	3.3 x 10 ⁵	-0.84	Acceptable Recovery
		2	2.3 x 10 ⁶	9.3 x 10 ⁵	2.9 x 10 ⁵	-0.90	Acceptable Recovery
		3	2.2 x 10 ⁶	9.3 x 10 ⁵	2.5 x 10 ⁵	-0.94	Acceptable Recovery
<i>Clostridium sporogenes</i> ATCC 3584	diluted 1:10	1	6.5 x 10 ⁵	2.7 x 10 ⁵	1.1 x 10 ⁵	-0.77	Acceptable Recovery
		2	6.4 x 10 ⁵	2.6 x 10 ⁵	9.9 x 10 ⁴	-0.81	Acceptable Recovery
		3	6.4 x 10 ⁵	2.6 x 10 ⁵	1.0 x 10 ⁵	-0.81	Acceptable Recovery
<i>Fusobacterium necrophorum</i> ATCC 25286	diluted 1:10	1	9.6 x 10 ⁵	2.7 x 10 ⁵	1.3 x 10 ⁵	-0.87	Acceptable Recovery
		2	9.7 x 10 ⁵	2.6 x 10 ⁵	1.2 x 10 ⁵	-0.91	Acceptable Recovery
		3	9.4 x 10 ⁵	2.6 x 10 ⁵	1.4 x 10 ⁵	-0.83	Acceptable Recovery
<i>Peptococcus magnus</i> ATCC 29328	diluted 1:10	1	4.9 x 10 ⁶	2.8 x 10 ⁶	6.9 x 10 ⁵	-0.85	Acceptable Recovery
		2	4.9 x 10 ⁶	2.7 x 10 ⁶	5.3 x 10 ⁵	-0.97	Acceptable Recovery
		3	4.8 x 10 ⁶	2.6 x 10 ⁶	5.7 x 10 ⁵	-0.93	Acceptable Recovery

In accordance with CLSI M40-A, with the exception of *Neisseria gonorrhoeae*, viability performance is measured for each test organism at the 48 h time point and compared with the acceptance criteria. Viability performance is measured for *Neisseria gonorrhoeae* at the 24 h time point. In both the Roll-Plate and Swab Elution viability performance studies, BD ESwab System was able to maintain acceptable recovery of all organisms evaluated at both refrigerator (4 – 8°C) and room temperature (20 – 25°C). Acceptable recovery for the Roll-Plate Method is defined as ≥ 5 CFU following the specified holding time from the specific dilution that yielded zero-time plate counts closest to 300 CFU. Acceptable recovery for the Swab Elution Method is defined as no more than a $3 \log_{10}$ ($1 \times 10^3 \pm 10\%$) decline in CFU between the zero-time CFU count and the CFU of the swabs after the specified holding time.

Viability performance studies also include an assessment of bacterial overgrowth at refrigerated temperatures (4 – 8°C). For the Swab Elution Method, an overgrowth assessment is made on all bacteria species tested at the 48 h holding time point except for *Neisseria gonorrhoeae* which is assessed at the 24 h holding time point. Overgrowth assessment using the Swab Elution Method is defined as greater than $1 \log_{10}$ increase in CFU between the zero-time CFU count and the holding time point. For the Roll-Plate Method, an overgrowth assessment is made with a separate analysis in which swabs are dosed with 100 μ L containing 10^2 CFU of *Pseudomonas aeruginosa* culture. Overgrowth under these conditions is defined as greater than $1 \log_{10}$ increase in CFU between zero-time CFU and the 48 h holding time point. BD ESwab Collection and Transport System demonstrated no overgrowth in either the Swab Elution or Roll-Plate Methods based on the acceptance criteria described in CLSI M40-A.

AVAILABILITY

Cat. No.	Description
220245	BD™ ESwab polypropylene screw-cap tube with 1 mL of Liquid Amies Medium with a white cap and one regular size applicator swab with flocked nylon fiber tip, 50 kits per shelf pack
220246	BD™ ESwab polypropylene screw-cap tube with 1 mL of Liquid Amies Medium with a green cap and one minitip swab with flocked nylon fiber tip, 50 kits per shelf pack
220532	BD™ ESwab polypropylene screw-cap tube with 1 mL of Liquid Amies Medium with a blue cap and one flexible minitip swab with flocked nylon fiber tip, 50 kits per shelf pack

REFERENCES

- Amies CR. 1967. A modified formula for the preparation of Stuart's medium. Canadian Journal of Public Health, Vol. 58, 296 – 300.
- Miller JM. 1999. A guide to specimen management in clinical microbiology, 2nd ed. American Society for Microbiology. Washington, DC.
- Miller JM, Holmes HT. 1995. Specimen collection, transport, and storage. In: Manual of clinical microbiology, 6th ed. (Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH, eds.), ASM, Washington, DC.
- Clinical and Laboratory Standards Institute. 2003. Quality control of microbiological transport systems; Approved Standard. M40-A Vol. 23 No. 34.
- Perry JL, Matthews JS. 2003. Compliance of two popular swab transport systems with performance standards detailed by the new NCCLS Proposed Standard, M40-P. 103rd General Meeting of the American Society for Microbiology, Washington, DC. Abstract C-42.
- Robinson A, Gruver ML. 2002. Comparison of bacterial survival in two transport systems stored at room temperature and refrigerator temperatures. 102nd General Meeting of the American Society for Microbiology, Salt Lake City, UT. Abstract C-69.
- Human RP, Jones GA. 2004. Evaluation of 4 transport systems against a published standard. 104th General Meeting of the American Society for Microbiology, New Orleans, LA. Abstract C-161.
- Human RP, Jones GA. 2004. Evaluation of swab transport systems against a published standard. J. Clin. Pathol. 57:762-763.
- Arbique J, Campbell S, MacFarlane M, Davidson RJ. 2003. Comparison of methodologies for anaerobic organisms described in NCCLS document M40-P. Quality Control of Microbiology Transport Devices. 13th European Congress of Clinical Microbiology and Infectious Disease (ECCMID), Glasgow, UK. Abstract P-652.
- Clinical and Laboratory Standards Institute. 2005. Approved Guideline M29-A3. Protection of laboratory workers from occupationally acquired infections, 3rd ed. CLSI, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17: 53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Koneman EW, Allen SD, Janda WM, Schreckenberger PC and Winn, Jr. WC. 1992. Color atlas and textbook of diagnostic microbiology. 4th ed. J.B. Lippincott Co. Philadelphia, PA.
- 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- Forbes BA, Sahm DF, Weissfeld AS. 1998. Bailey and Scott's Diagnostic Microbiology. 10th ed. Mosby, St. Louis, MO.
- Isenberg HD. 2004. Clinical microbiology procedures handbook, 2nd ed. ASM, Washington, DC.
- Isenberg HD. 1998. Essential procedures for clinical microbiology. Chapter 14.12, Page 787. Packaging and shipping infectious substances. ASM, Washington, DC.
- Clinical and Laboratory Standards Institute. 1994. Procedures for handling and transport of diagnostic specimens and etiologic agents; Approved Standard. H5-A3.
- Isenberg HD, Schoenkench FD, Von Graevenitz A. 1979. Cumitech 9, Collection and processing of bacteriological specimens. Coordinating editor, SJ. Rubin. American Society for Microbiology, Washington, DC.

21. Murray PR, Baron EJ, Pfaller MA, Tenover FC, Tenover RH. 1999. Manual of clinical microbiology. 7th ed., ASM Washington, DC.
22. Summanen P, Baron EJ, Citron D, Strong C, Wexler HM, Finegold SM. 1993. Wadsworth anaerobic bacteriology manual, 5th ed. Star Publishing Company, Belmont, CA.
23. Isenberg HD. 1998. Essential procedures for clinical microbiology. Chapter 1.1, Page 27. Collection, transport and manipulation of clinical specimens. Procedure for streaking plates for primary isolation. ASM, Washington, DC.
24. Marler LM, Siders JA, Allen SD. 2001. Direct smear atlas, a monograph of gram-stained preparations of clinical specimens. Lippincott Williams and Wilkins.
25. Rotimi VO, Yakubu Z, Abudu OO, Banjo TO. 1991. Direct Gram's stain of vaginal discharge as a means of diagnosing bacterial vaginosis. Journal of Medical Microbiology, 35:103-106.
26. Spiegel CA, Amsel R, Holmes KK. 1983. Diagnosis of bacterial vaginosis by direct Gram stain of vaginal fluid J. Clin Microbiol. 18:170-177.
27. Benavides MI, Moncada X, Rodriguez B, Castillo C. 1992. Gonococcal urethritis in men: clinical experience in 1978-1988. Rev Med Chil. 120(10):1140-3.
28. Mayaud P, Msuya W, Todd J, Kaatano G, West B, Begkoyan G, Grosskurth H, Mabey D. 1997. Rapid assessment in Rwandan refugee camps in Tanzania. Genitourin Med, Feb;73 (1):33-8.
29. Deceuninck G, Asamoah-Adu C, Khonde N, Pepin J, Frost EH, Deslandes S, Asamoah-Abu A, Bekoe V, Alary M. 2000. Improvement of clinical algorithms for the diagnosis of *Neisseria gonorrhoeae* and *Chlamydia trachomatis* by the use of Gram-stained smears among female sex workers in Accra, Ghana. Sex Transm Dis. Aug;27 (7):401-10.
30. Washington JA. 1986. Rapid diagnosis by microscopy. Clin. Microbiol. News. 8:135-137.
31. Isenberg HD. 1998. Essential procedures for clinical microbiology. Chapter 2.1, Page 41. Gram stain. ASM, Washington, DC.
32. Arbiq J, Campbell S, MacFarlane M, Davidson RJ. 2003. Comparison of methodologies described in NCCLS document M40-P Quality Control of Microbiology Transport Devices. 103rd General Meeting of the American Society for Microbiology, Washington, DC. Abstract C-40.
33. Van Horn KG, Rankin I. 2005. Evaluation and comparison of two Stuart's Liquid Swab transport systems tested by the NCCLS M40 method. 105th General Meeting of the American Society for Microbiology, Atlanta, Georgia. Abstract C-292.
34. Bourbeau PP, Heiter BJ. 2003. Validation of QC standard for bacteriological transport devices as specified in the NCCLS Proposed Standard M40: Quality Control of Microbiological Transport Systems. 103rd General Meeting of the American Society for Microbiology, Washington, DC. Abstract C-46.

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Becton Dickinson France S.A.S.